

m INIPHOR TS

Electrophoresis Automated System
High Technology in Italian Style Instrument



SELEO SRL



Meridian Healthcare®
medical diagnostics and equipment

mINIPHOR TS

Electrophoresis Automated System

High Technology in Italian Style Instrument

MINIPHOR TS is the latest of the family of electrophoresis instruments; the smallest fully automatic instrument, with an integrated control system and a touch-screen interface to view all the phases of process execution and the correction and memorization of pherograms.

MINIPHOR TS allows to perform the most widespread electrophoretic techniques: Serumproteins, Haemoglobins, Lipoproteins, Urinary Proteins.

Touch-screen interface with Windows IoT Enterprise operating system, LAN connection for remote host protocol.

Automatic handling of reagents with washing procedure.

YOUR PARTNER IN ELECTROPHORESIS



Cellulose Acetate membranes
Balanced kits and applicative support



DISTRIBUITO IN ITALIA DA

MERIDIAN HEALTHCARE srl

Via Caronda, 446 SC/A - 95100 CATANIA - ITALIA

Tel.: +39 095 7256869 Fax: +39 095 7254454

e-mail: info@meridianhealthcare.it

web: www.meridianhealthcare.it