


EmA ENDOMYSIUM IgG ANTIBODIES

EmA human ELISA kit is a quantitative immunoassay for the determination of IgG antibodies against endomysial antigens (EmA) in human serum. EmA human ELISA kit is intended as an aid in the diagnosis of diagnosis of celiac disease and dermatitis herpetiformis in conjunction with other clinical and laboratory findings.



- Standardized ELISA test format.
 - Ready to use reagents.
- Microtiter plate coated with with purified human endomysial antigens and gliadin fragments.
 - 🕒 Incubation: 105 min.
- Positive Control: ready to use
- Negative Control: ready to use
 - Sensitivity: 95.2%
 - Specificity: 86.3%
 - Test: 96

ORDERING INFORMATION:



Code No.	Product	Pack-size
4045	Anti-Endomysium IgG	96 wells

EmA anti-ENDOMISIO IgG ELISA

Il kit ELISA EmA è un test immunologico quantitativo per la determinazione degli anticorpi IgG contro gli antigeni endomysiali (EmA) nel siero umano. Il dispositivo è inteso come ausilio nella diagnosi della malattia celiaca e della dermatite erpetiforme insieme ad altri risultati clinici e di laboratorio.

- Kit ELISA standardizzato.
- Reagenti liquidi e pronto all'uso.
- Piastra per microtitolazione rivestita con antigeni endomysiali umani purificati e frammenti di gliadina.
 - 🕒 Incubazione: 105 min.
- Controllo Positivo: pronto all'uso.
- Controllo Negativo: pronto all'uso.
 - Sensibilità: 95.2%
 - Specificità: 86.3%
 - Test: 96

